

CHAMBER PLUG CONFIGURATIONS:

MODEL	DESCRIPTION	VOLTAGE	PLUG TYPE
2050-1-1	Refrigerated & Heated Bath/Circulator	115V, 60 Hz	NEMA 5-15
2050-1-4	Refrigerated & Heated Bath/Circulator	230V, 50/60 Hz	CEE 7/7
6105-2	Fingerprint Development Chamber, 5 cu. ft.	208/230V, 60 Hz	NEMA 6-15
6400-1	Oasis CO2 Incubator, 6 cu. ft.	115V, 60 Hz	NEMA 5-15
6400-4	Oasis CO2 Incubator, 6 cu. ft.	230V, 50/60 Hz	CEE 7/7
6404-1	Oasis CO2/O2 Incubator, 6 cu. ft.	115V, 60 Hz	NEMA 5-15
6404-4	Oasis CO2/O2 Incubator, 6 cu. ft.	230V, 50/60 Hz	CEE 7/7
CRSY102-1	Universal Condensate Recirculator	115V, 60 Hz	NEMA 5-15
CRSY102-4	Universal Condensate Recirculator	230V, 50/60 Hz	CEE 7/7
GASG201-1	Stand Alone Gas Guard/Tank Switcher	115V, 60 Hz	NEMA 5-15
GASG201-4	Stand Alone Gas Guard/Tank Switcher	230V, 50/60 Hz	NEMA 6-15
7000-10-1	Environmental Chamber, 10 cu. ft.	115V, 60 Hz	NEMA 5-15
7000-10-2	Environmental Chamber, 10 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7000-25-1	Environmental Chamber, 25 cu. ft.	115V, 60 Hz	NEMA 5-20
7000-25-2	Environmental Chamber, 25 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7000-33-1	Environmental Chamber, 33 cu. ft.	115V, 60 Hz	NEMA 5-20
7000-33-2	Environmental Chamber, 33 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7000-50-1	Environmental Chamber, 2-Door, 50 cu. ft.	115V, 60 Hz	NEMA 5-20
7000-50-2	Environmental Chamber, 2-Door, 50 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7000-75-1	Environmental Chamber, 3-Door, 75 cu. ft.	115V, 60 Hz	NEMA 5-20
7000-75-2	Environmental Chamber, 3-Door, 75 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7001-10-1	Refrigerated Incubator, 10 cu. ft.	115V, 60 Hz	NEMA 5-15
7001-10-2	Refrigerated Incubator, 10 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7001-25-1	Refrigerated Incubator, 25 cu. ft.	115V, 60 Hz	NEMA 5-20
7001-25-2	Refrigerated Incubator, 25 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7001-33-1	Refrigerated Incubator, 33 cu. ft.	115V, 60 Hz	NEMA 5-20
7001-33-2	Refrigerated Incubator, 33 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7001-50-1	Refrigerated Incubator, 2-Door, 50 cu. ft.	115V, 60 Hz	NEMA 5-20
7001-50-2	Refrigerated Incubator, 2-Door, 50 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7001-75-1	Refrigerated Incubator, 3-Door, 75 cu. ft.	115V, 60 Hz	NEMA 5-20
7001-75-2	Refrigerated Incubator, 3-Door, 75 cu. ft.	208/230V, 60 Hz	NEMA 6-15

MODEL	DESCRIPTION	VOLTAGE	Plug Type
7002-10-1	Heated Incubator, 10 cu. ft.	115V, 60 Hz	NEMA 5-15
7002-10-2	Heated Incubator, 10 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7002-25-1	Heated Incubator, 25 cu. ft.	115V, 60 Hz	NEMA 5-20
7002-25-2	Heated Incubator, 25 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7002-33-1	Heated Incubator, 33 cu. ft.	115V, 60 Hz	NEMA 5-20
7002-33-2	Heated Incubator, 33 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7003-10-1	Humidified Incubator, 10 cu. ft.	115V, 60 Hz	NEMA 5-15
7003-10-2	Humidified Incubator, 10 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7003-25-1	Humidified Incubator, 25 cu. ft.	115V, 60 Hz	NEMA 5-20
7003-25-2	Humidified Incubator, 25 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7003-33-1	Humidified Incubator, 33 cu. ft.	115V, 60 Hz	NEMA 5-20
7003-33-2	Humidified Incubator, 33 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7300-22-1	Plant Growth, 1-Door, 22 cu. ft.	115V, 60 Hz	NEMA 5-20
7300-22-2	Plant Growth, 1-Door, 22 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7300-50-2	Plant Growth, 2-Door, 50 cu. ft.	208/230V, 60 Hz	NEMA 6-20
7300-75-2	Plant Growth, 3-Door, 75 cu. ft.	208/230V, 60 Hz	NEMA L6-30
7301-22-1	Plant Growth, 1-Door, 22 cu. ft.	115V, 60 Hz	NEMA 5-20
7301-22-2	Plant Growth, 1-Door, 22 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7301-50-2	Plant Growth, 2-Door, 50 cu. ft.	208/230V, 60 Hz	NEMA 6-20
7301-75-2	Plant Growth, 2-Door, 50 cu. ft.	208/230V, 60 Hz	NEMA L6-30
7303-22-1	Arabidopsis/Algae, 1-Door, 22 cu. ft.	115V, 60 Hz	NEMA 5-20
7303-22-2	Arabidopsis/Algae, 1-Door, 22 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7303-50-2	Arabidopsis/Algae, 2-Door, 50 cu. ft.	208/230V, 60 Hz	NEMA 6-20
7303-75-2	Arabidopsis/Algae, 3-Door, 75 cu. ft.	208/230V, 60 Hz	NEMA L6-30
7304-22-1	Arabidopsis/Algae, 1-Door, 22 cu. ft.	115V, 60 Hz	NEMA 5-20
7304-22-2	Arabidopsis/Algae, 1-Door, 22 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7304-50-2	Arabidopsis/Algae, 2-Door, 50 cu. ft.	208/230V, 60 Hz	NEMA 6-20
7304-75-2	Arabidopsis/Algae, 3-Door, 75 cu. ft.	208/230V, 60 Hz	NEMA L6-30
7305-22-1	Plant Tissue Culture, 1-Door, 22 cu. ft.	115V, 60 Hz	NEMA 5-20
7305-22-2	Plant Tissue Culture, 1-Door, 22 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7306-22-1	Seed Germination, 1-Door, 22 cu. ft.	115V, 60 Hz	NEMA 5-20
7306-22-2	Seed Germination, 1-Door, 22 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7310-22-1	LED Plant Growth, 1-Door, 22 cu. ft.	115V, 60 Hz	NEMA 5-15
7310-22-2	LED Plant Growth, 1-Door, 22 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7310-50-2	LED Plant Growth, 2-Door, 50 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7310-75-2	LED Plant Growth, 3-Door, 75 cu. ft.	208/230V, 60 Hz	NEMA 6-15

MODEL	DESCRIPTION	VOLTAGE	Plug Type
7311-22-1	LED Plant Growth, 1-Door, 22 cu. ft.	115V, 60 Hz	NEMA 5-15
7311-22-2	LED Plant Growth, 1-Door, 22 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7311-50-2	LED Plant Growth, 2-Door, 50 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7311-75-2	LED Plant Growth, 3-Door, 75 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7312-22-1	LED Plant Growth, 1-Door, 22 cu. ft.	115V, 60 Hz	NEMA 5-15
7312-22-2	LED Plant Growth, 1-Door, 22 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7312-50-2	LED Plant Growth, 2-Door, 50 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7312-75-2	LED Plant Growth, 3-Door, 75 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7315-22-1	LED Plant Growth, 1-Door, 22 cu. ft.	115V, 60 Hz	NEMA 5-15
7315-22-2	LED Plant Growth, 1-Door, 22 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7317-22-1	LED Plant Growth, 1-Door, 22 cu. ft.	115V, 60 Hz	NEMA 5-15
7317-22-2	LED Plant Growth, 1-Door, 22 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7317-50-2	LED Plant Growth, 2-Door, 50 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7317-75-2	LED Plant Growth, 3-Door, 75 cu. Ft	208/230V, 60 Hz	NEMA 6-15
7340-25-1	Insect Rearing Chamber, 25 cu. ft.	115V, 60 Hz	NEMA 5-20
7340-25-2	Insect Rearing Chamber, 25 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7340-33-1	Insect Rearing Chamber, 33 cu. ft.	115V, 60 Hz	NEMA 5-20
7340-33-2	Insect Rearing Chamber, 33 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7350-25-1	Animal Rearing Chamber, 25 cu. ft.	115V, 60 Hz	NEMA 5-20
7350-25-2	Animal Rearing Chamber, 25 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7350-33-1	Animal Rearing Chamber, 33 cu. ft.	115V, 60 Hz	NEMA 5-20
7350-33-2	Animal Rearing Chamber, 33 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7400-25-1	CO2 Incubator, Reach-In, 25 cu. ft.	115V, 60 Hz	NEMA 5-20
7400-25-2	CO2 Incubator, Reach-In, 25 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7400-33-1	CO2 Incubator, Reach-In, 33 cu. ft.	115V, 60 Hz	NEMA 5-20
7400-33-2	CO2 Incubator, Reach-In, 33 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7401-25-1	CO2 Incubator, Shaker Ready, 25 cu. ft.	115V, 60 Hz	NEMA 5-20
7401-25-2	CO2 Incubator, Shaker Ready, 25 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7401-33-1	CO2 Incubator, Shaker Ready, 33 cu. ft.	115V, 60 Hz	NEMA 5-20
7401-33-2	CO2 Incubator, Shaker Ready, 33 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7402-25-1	CO2 Incubator, Cell Roller Ready, 25 cu. ft.	115V, 60 Hz	NEMA 5-20
7402-25-2	CO2 Incubator, Cell Roller Ready, 25 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7402-33-1	CO2 Incubator, Cell Roller Ready, 33 cu. ft.	115V, 60 Hz	NEMA 5-20
7402-33-2	CO2 Incubator, Cell Roller Ready, 33 cu. ft.	208/230V, 60 Hz	NEMA 6-15

MODEL	DESCRIPTION	VOLTAGE	Plug Type
7403-25-1	CO2 Incubator, Shaker/Cell Roller Ready, 25 cu. Ft.	115V, 60 Hz	NEMA 5-20
7403-25-2	CO2 Incubator, Shaker/Cell Roller Ready, 25 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7403-33-1	CO2 Incubator, Shaker/Cell Roller Ready, 33 cu. ft.	115V, 60 Hz	NEMA 5-20
7403-33-2	CO2 Incubator, Shaker/Cell Roller Ready, 33 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7404-10-1	CO2 Incubator, GelJacket, 10 cu. ft.	115V, 60 Hz	NEMA 5-15
7404-10-2	CO2 Incubator, GelJacket, 10 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7410-5-1	Wally GT CO2 Incubator, 5 cu. ft.	115V, 60 Hz	NEMA 5-15
7410-5-2	Wally GT CO2 Incubator, 5 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7411-5-1	Wally GT CO2/O2 Incubator, 5 cu. ft.	115V, 60 Hz	NEMA 5-15
7411-5-2	Wally GT CO2/O2 Incubator, 5 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7420-5-1	Wally CO2 Incubator, 5 cu. ft.	115V, 60 Hz	NEMA 5-15
7420-5-2	Wally CO2 Incubator, 5 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7540-11-1	Photostability Chamber, 11 cu. ft.	115V, 60 Hz	NEMA 5-20
7540-11-2	Photostability Chamber, 11 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7545-11-1	Photostability Chamber, 11 cu. ft.	115V, 60 Hz	NEMA 5-20
7545-11-2	Photostability Chamber, 11 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7900-25-2	Freeze/Thaw Chamber, 25 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7900-33-2	Freeze/Thaw Chamber, 33 cu. Ft.	208/230V, 60 Hz	NEMA 6-15
7901-25-2	Freeze/Thaw Chamber, 25 cu. ft.	208/230V, 60 Hz	NEMA 6-15
7901-33-2	Freeze/Thaw Chamber, 33 cu. ft.	208/230V, 60 Hz	NEMA 6-15

PLUG PICTURES

