

SPECIFICATIONS

MODEL	7000-50-1	7000-50-2	7000-50-3	7000-75-1	7000-75-2	7000-75-3
Temperature Range	5°C to 70°C					
Temperature Control	± 0.2°C					
Temperature Uniformity	± 1.0°C					
Temperature Sensor	3-wire RTD					
Humidity Range	20-95% RH					
Humidity Control	± 2%					
Humidity Sensor	Capacitive					
Interior Dimensions	51" W x 28.5" D x 59" H (130 cm x 72 cm x 150 cm)			79.5" W x 28.5" D x 59" H (202 cm x 72 cm x 150 cm)		
Interior Construction	Type 304, Stainless Steel					
Exterior Dimensions	57" W x 35.1" D* x 81.5" H (145 cm x 89 cm x 207 cm)			85.5" W x 35.1" D* x 81.5" H (217 cm x 89 cm x 207 cm)		
Exterior Construction	Stainless Steel					
Work Space	50 cu. ft. (1,416 Liters)			75 cu. ft. (2,124 Liters)		
# of Shelves	8			12		
Shelf Construction	Type 304, Perforated Stainless Steel, Electropolished					
Shelf Dimensions	21.6" W x 26.5" D (55 cm x 67 cm)					

Model	7000-50			7000-75		
	-1	-2	-3	-1	-2	-3
Shipping Weight	800 lbs (363 kg.)		1,100 lbs. (499 kg.)**	1125 lbs. (510 kg.)		1650 lbs. (748 kg.)**
Electrical	115V, 60Hz, 16A	230V, 60Hz, 12A	230V, 50Hz, 10A	115V, 60Hz, 16A	230V, 60Hz, 12A	230V, 50Hz 10A

Specifications are subject to change without notice.

Environmental Conditions: Temperature 15°C to 25°C, Humidity non-condensing

*See graph for details

**Includes export shipping crate