

# ANIMAL REARING CHAMBERS

*Models: 7350-25 / 7350-33*

Model	7350-25-1	7350-25-2	7350-25-3	7350-33-1	7350-33-2	7350-33-3
Temperature Range	10°C to 40°C					
Temperature Control	±0.1°C					
Temperature Uniformity	±0.3°C					
Lights	(2) Side-Mounted LED Modules					
Interior Dimensions	32" W x 27" D x 52.7"H (81.3 cm x 68.6 cm x 133.9 cm)			32" W x 27" D x 65.7"H (81.3 cm x 68.6 cm x 166.9 cm)		
Interior Construction	2B Finish Solid Stainless Steel (Type 304)					
Exterior Dimensions	35.5" W x 33.3" D* x 77.1" H (90.2 cm x 84.6 cm x 195.8 cm)			35.5" W x 33.3" D* x 90.1" H (90.2 cm x 84.6 cm x 228.9 cm)		
Exterior Construction	Powder Coated Cold Rolled Steel					
Work Space	25 ft <sup>3</sup> (708 liters)			33 ft <sup>3</sup> (934 liters)		
Shelf Count	4 Standard / 25 Maximum			5 Standard / 31 Maximum		
Shelf Construction	Perforated Stainless Steel (Type 304)					
Shelf Dimensions	29.25" W x 24.45" D (74.3 cm x 62.1 cm)			29.25" W x 24.45" D (74.3 cm x 62.1 cm)		
Electrical	115V 60 Hz 16A	208/230V 60Hz 10A	230V 50 Hz 8A	115V 60 Hz 16A	208/230V 60 Hz 10A	230V 50 Hz 8A
Shipping Weight lbs.	725 lbs.		825 lbs.**	775 lbs.		1,050 lbs.**
Shipping Weight kg.	329 kg.		374 kg.**	352 kg.		476 kg.**

Specifications are based on 20°C ambient and standard voltage. Specifications are subject to change without notice.

\*Add 2.75" for handle. \*\*7350-25-3 and 7350-33-3; Includes Export Shipping Crate.