

CO₂ INCUBATOR / GELJACKET

Model: 7404-10

Model	7404-10-1	7404-10-2	7404-10-3
Temperature Range	10°C above ambient to 60°C		
Temperature Control	±0.1°C		
Temperature Uniformity	±0.3°C		
Temperature Sensor	RTD		
Humidity Range	Ambient to 95%RH		
Humidity Control	±3% RH		
Humidity Sensor	Capacitive		
CO ₂ Range	0-20% CO ₂		
CO ₂ Control	±0.1% CO ₂		
CO ₂ Sensor	Infrared CO ₂ Sensor		
Interior Dimensions	23" W x 25.8" D x 29.8" H (58.4 cm x 65.5 cm x 75.7 cm)		
Interior Construction	Type 304, 2B Finish, Solid Stainless Steel		
Exterior Dimensions	44.5" W x 31.7" D* x 37" H (112.3 cm x 82.8 cm x 92.7 cm)		
Exterior Construction	Cold Rolled Steel, Powder Coated		
Work Space	10 cu. ft. (283 Liters)		
# of Shelves	3 Standard; 15 Maximum		
Shelf Construction	Type 304, Perforated Stainless Steel, Electropolished		
Shelf Dimensions	22.5" W x 25.4" D (57.2 cm x 64.5 cm)		
Electrical	115V, 60 Hz, 12A	208/230V, 60 Hz, 8A	230V, 50 Hz, 6A
Shipping Weight	450 lbs.		750 lbs.**

Specifications are subject to change without notice. *Add 2.75" for handle. ** Includes export shipping crate.