

CO₂ INCUBATOR / REACH-IN

Models: 7400-25 / 7400-33

Model	7400-25-1	7400-25-2	7400-25-3	7400-33-1	7400-33-2	7400-33-3
Temperature Range	10°C above ambient to 60°C					
Temperature Control	±0.1°C					
Temperature Uniformity	±0.3°C					
Temperature Sensor	RTD					
Humidity Range (with standard water pan)	Elevated up to 95% @ 37°C					
Controlled Humidity Range (with optional HUMD304)	Ambient to 95% RH, ±3% RH					
CO ₂ Range	0-20% CO ₂					
CO ₂ Control	±0.1% CO ₂					
CO ₂ Sensor	Infrared CO ₂ Sensor					
Interior Dimensions	32" W x 27" D x 52.7" H (81.3 cm x 68.6 cm x 133.9 cm)			32" W x 27" D x 65.7" H (81.3 cm x 68.6 cm x 166.9 cm)		
Interior Construction	Type 304, 2B Finish Solid Stainless Steel					
Exterior Dimensions	35.5" W x 33.3" D* x 77.1" H (90.2 cm x 84.6 cm x 195.8 cm)			35.5" W x 33.3" D* x 90.1" H (90.2 cm x 84.6 cm x 228.9 cm)		
Exterior Construction	Powder Coated Cold Rolled Steel					
Work Space	25 ft ³ (708 liters)			33 ft ³ (934 Liters)		
# of Shelves	4 Standard; 25 Maximum			5 Standard; 31 Maximum		
Shelf Construction	Type 304, Perforated Stainless Steel, Electropolished					
Shelf Dimensions	29.25" W x 24.45" D (74.3 cm x 62.1 cm)					
Electrical	115V 60 Hz 16A	208/230V 60 Hz 10A	230V 50 Hz 10A	115V 60 Hz 16A	208/230V 60 Hz 10A	230V 50 Hz 10A
Shipping Weight lbs.	725		825**	775		1,050**
Shipping Weight kg.	329		374**	352		476**

Specifications are based on 20°C ambient and standard voltage. Specifications are subject to change without notice.

*Add 2.75" for handle. **Includes Export Shipping Crate.